



3Dconnexion Devices Certification *

Release 2021 R1

* Certification means that the 3Dconnexion products listed work as expected with the current release of the ANSYS applications listed below.

- Use of the most recent driver version is recommended for the listed operating systems.
- For help with ANSYS products, please contact your ANSYS support provider.
- For device support, see www.3dconnexion.com.

| | | 3Dconnexion Product | | | | Tested Operating Systems |
|--------------------------------------------------------------------------------------------|---------------------------------|---------------------|---------------------|------------------|-----------------------|--------------------------|
| | | SpaceNavigator | SpaceMouse Wireless | SpaceMouse Pro | SpaceMouse Enterprise | |
| Workbench-hosted Applications for interactive desktop execution | ANSYS CFX | ✓ ¹ | ✓ ¹ | ✓ ¹ | ✓ ¹ | Windows 10 |
| | ANSYS Discovery | ✓ | ✓ | ✓ | ✓ | Windows 10 |
| | ANSYS Discovery SpaceClaim | ✓ | ✓ | ✓ | ✓ | Windows 10 |
| | ANSYS ICEM CFD | ✓ ¹ | ✓ ¹ | ✓ ^{1,4} | ✓ ^{1,4} | Windows 10 |
| Stand-alone Applications for interactive or batch execution outside of Workbench | ANSYS Icepak | ✓ ^{1,3} | ✓ ^{1,3} | ✓ ^{1,3} | ✓ ^{1,3} | Windows 10 |
| | ANSYS Mechanical | ✓ | ✓ | ✓ ² | ✓ | Windows 10 |
| | ANSYS Mechanical APDL | ✓ | ✓ | ✓ | ✓ | Windows 10 |
| | ANSYS SpaceClaim Direct Modeler | ✓ | ✓ | ✓ | ✓ | Windows 10 |

Product Functionality

1. View buttons are not functional.
2. Lock Rotation button not functional.
3. Redo button not functional.
4. Rotate button works as fit.